

RPCB1493

Product Information

Product SKU:	RPCB1493	Gene ID:	1490	Size:	10µg
Tag:	No tag	Reactivity:	Human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	P29279
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: This protein is a mitogen that is secreted by vascular endothelial cells. The encoded protein plays a role in chondrocyte proliferation and differentiation, cell adhesion in many cell types, and is related to platelet-derived growth factor. Certain polymorphisms in this gene have been linked with a higher incidence of systemic sclerosis.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IGFBP-8/CTGF Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH7.4.Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.