

RPCB1515

---

**Product Information**

<b>Product SKU:</b>	RPCB1515	<b>Gene ID:</b>	43740568	<b>Size:</b>	100µg
<b>Tag:</b>	C-His&Avi	<b>Reactivity:</b>	SARS-CoV-2		

---

**Additional Information**

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P0DTC2
<b>Purity:</b>	> 95% by SDS-PAGE.> 95% by HPLC.		

---

**Protein Information**

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV-2 Spike RBD Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1.0 EU/µg of the protein by LAL method.

**Formulation:** Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4 or Supplied as a 0.22 µm filtered solution in PBS, pH 7.4.Contact us for customized product form or formulation.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.