

# Recombinant SARS-CoV-2 Spike S1 (D614G) Protein RPCB1516



## Product Information

<b>Product SKU:</b>	RPCB1516	<b>Gene ID:</b>	43740568	<b>Size:</b>	50µg
<b>Tag:</b>	C-His	<b>Reactivity:</b>	SARS-CoV-2		

## Additional Information

<b>Expression Host:</b>	HEK293 cells	<b>Swissprot:</b>	P0DTC2
<b>Purity:</b>	> 95% by SDS-PAGE.		

## Protein Information

**Background:** -

**Protein Description:** High quality, high purity and low endotoxin recombinant Recombinant SARS-CoV-2 Spike S1(D614G) Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

**Endotoxin:** < 1.0 EU/µg of the protein by LAL method.

**Formulation:** Supplied as a 0.22 µm filtered solution in PBS, pH 7.4 or Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

**Storage:** Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.

**Contact Details | Dublin, Ireland**

**Email:** techsupport@assaygenie.com | **Web:** [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.