

RPCB1520

Product Information

Product SKU:	RPCB1520	Gene ID:	146433	Size:	10µg
Tag:	C-6His	Reactivity:	human		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	Q6ZMJ4
Purity:	> 95% by SDS-PAGE.		

Protein Information

Background: Interleukin are a group of cytokines that play an important role in the immune system. They modulate inflammation and immunity by regulating growth, mobility and differentiation of lymphoid and other cells. This entry represents interleukin-34 (IL-34), it was identified via functional screening of a library of secreted proteins. This cytokine promotes the differentiation and viability of monocytes and macrophages through the colony-stimulating factor-1 receptor (CSF1R)

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human IL-34 Protein, tested reactivity in HEK293 cells and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: < 0.01EU/µg

Formulation: Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.

Storage: Store at -20°C. Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt. After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.