

RPCB1603

Product Information

Product SKU:	RPCB1603	Gene ID:	27252	Size:	100µg
Tag:	No-tag	Reactivity:	Human		

Additional Information

Expression Host:	-	Swissprot:	Q9Y2M5
Purity:	-		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human KLHL20 Protein , tested reactivity inE. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: Please contact us for more information.

Formulation: Supplied as a 0.22 µm filtered solution in PBS pH7.4.

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.