

RPCB1742

Product Information

Product SKU:	RPCB1742	Gene ID:	6281	Size:	10µg
Tag:	C-6His	Reactivity:	Human		

Additional Information

Expression Host:	E. coli	Swissprot:	P60903
Purity:	-		

Protein Information

Background: S100A10 functions as a linker tethering certain transmembrane proteins to annexin A2 thereby assisting their traffic to the plasma membrane and/or their firm anchorage at certain membrane sites.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human S100-A10 Protein , tested reactivity in E. coli and has been validated in SDS-PAGE. 100% guaranteed.

Endotoxin: <0.1 EU/µg

Formulation: Supplied as a 0.22 µm filtered solution in 20mM Tris 10% Glycerol, PH8.5

Storage: Store at -70°C. This product is stable at ≤ -70°C for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles.