

RPCB2002

Product Information

Product SKU:	RPCB2002	Gene ID:	57486	Size:	50µg
Tag:	N-His	Reactivity:	Human		

Additional Information

Expression Host:	E.coli	Swissprot:	Q9BYT8
Purity:	> 90% by SDS-PAGE.		

Protein Information

Background: -

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Human Neurolysin/NLN Protein , tested reactivity in E.coli and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: < 1.0 EU/µg of the protein by LAL method.

Formulation: Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization. Contact us for customized product form or formulation.

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.