

RPCB2091

Product Information

Product SKU:	RPCB2091	Gene ID:	102119254	Size:	100µg
Tag:	C-His	Reactivity:	Cynomolgus		

Additional Information

Expression Host:	HEK293 cells	Swissprot:	A0A2K5TXV0-1
Purity:	> 95% as determined by Tris-Bis PAGE ; > 95% as determined by HPLC		

Protein Information

Background: Thymic Stromal Lymphopoietin (TSLP) was originally identified as an activity from the conditioned medium of a mouse thymic stromal cell line that promoted the development of B cells. Cytokine that induces the release of T-cell-attracting chemokines from monocytes and, in particular, enhances the maturation of CD11c dendritic cells. Can induce allergic inflammation by directly activating mast cells.

Protein Description: High quality, high purity and low endotoxin recombinant Recombinant Cynomolgus TSLP Protein , tested reactivity in HEK293 cells and has been validated in SDS-PAGE.100% guaranteed.

Endotoxin: Less than 1EU per µg by the LAL method.

Formulation: -

Storage: Store at -20°C.Store the lyophilized protein at -20°C to -80 °C up to 1 year from the date of receipt.After reconstitution, the protein solution is stable at -20°C for 3 months, at 2-8°C for up to 1 week.