

# Recombinant SARS-CoV-2 Papain-Like Protease Protein

RPES0008

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES0008

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** SARS-CoV-2

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** PL-PRO, PLpro, Papain-like Protease, Replicase polyprotein 1a, pp1a

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** None

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 35.8 kDa

**Formulation:** Supplied as a 0.2 µM filtered solution of 20mM Tris-HCl, 10 mM 2-Mercaptoethanol, 20% Glycerol, pH 7.5.

**Observed MW:** 34 kDa

**Reconstitution:** -

**Accession:** QHD43415.1

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived SARS-CoV-2 Papain-Like Protease protein Glu1564-Lys1878

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Glu1564-Lys1878

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid