

Recombinant Human ACE2 Protein (Avi-His Tag)

RPES0029

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0029

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 100µg, 20µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: ACE-Related Carboxypeptidase, ACE2, ACEH, Angiotensin-Converting Enzyme 2, Angiotensin-Converting Enzyme Homolog, Metalloprotease MPROT15

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-Avi-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Immobilized 2019-nCoV S1 Protein-Fc(Cat#RPES0016) at 2µg/ml (100 µl/well) can bind Biotinylated Human ACE-2-Avi-His (Cat#RPES0029). The ED50 of Biotinylated Human ACE-2-Avi-His (Cat#RPES0029) is 0.17 ug/ml.

Calculated MW: 86.2 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 300mM NaCl, 1mM ZnCl₂ , 10% Glycerol, pH 7.4.

Observed MW: 110-140 kDa

Reconstitution: -

Accession: Q9BYF1

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Source: HEK293 Cells-derived Human ACE2 protein Gln18-Ser740, with an C-terminal Avi & His

Sequence: Gln18-Ser740

Form: Liquid

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.