

Recombinant Human NGAL/Lipocalin-2 Protein (His Tag, Human Cells)

RPES0435

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0435

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: 25 kDa alpha-2-microglobulin-related subunit of MMP-9, HNL, Lipocalin-2, NGAL, Neutrophil gelatinase-associated lipocalin, Oncogene 24p3, Siderocalin LCN2, p25

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 21.6 kDa

Formulation: Supplied as a 0.2 µm filtered solution of PBS, 50% Glycerol, pH 7.4.

Observed MW: 23 kDa

Reconstitution: -

Accession: P80188

Source: HEK293 Cells-derived Human NGAL, Lipocalin-2 protein Gln21-Gly198, with an C- terminal His

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Sequence: Gln21-Gly198

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid