

Recombinant Human HPGDS/GSTS Protein

RPES0453

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0453

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: GST Class-Sigma, GSTS, Glutathione S-Transferase, Glutathione-Dependent PGD Synthase, Glutathione-Requiring Prostaglandin D Synthase, H-PGDS, HPGDS, Hematopoietic Prostaglandin D Synthase, PGDS, PTGDS2, Prostaglandin-H2 D-Isomerase

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 22.3 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 200mM NaCl, pH 7.0.

Observed MW: 26 kDa

Reconstitution: -

Accession: O60760

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human HPGDS, GSTS protein Met 1-Leu199

Sequence: Met 1-Leu199

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid