

Recombinant Mouse EphB1/EPHT2 Protein (His & GST Tag)

RPES0628

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0628

Purity: > 80 % as determined by reducing SDS-PAGE.

Contents: 50µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: 9330129L11, AW488255, C130099E04Rik, Cek6, ENSMUSG00000074119, Elk, Elkh, Hek6, Net

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His-GST

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: Baculovirus-Insect Cells

Bio-Activity: The specific activity was determined to be 234 nmol/min/mg using Poly(Glu, Tyr) 4:1 as substrate.

Calculated MW: 72.4 kDa

Formulation: Supplied as sterile solution of 20mM Tris, 500mM NaCl, pH 8.0, 10% glycerol

Observed MW: 60 kDa

Reconstitution: -

Accession: Q8CBF3-1

Source: Baculovirus-Insect Cells-derived Mouse EphB1/EPHT2 protein Met 591-Ala 984, with an N-terminal His & GST

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Met 591-Ala 984

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid