

# Recombinant Mouse CD32b/FCGR2B Protein (His Tag)

RPES0629

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES0629

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Mouse

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** CD32, Fc gamma receptor IIB, Fc-gamma RII, Fc-gamma-RIIB, FcRII, Fcgr2, Fcgr2b, IgG Fc receptor II beta, Low affinity immunoglobulin gamma Fc region receptor II, Ly-17, Lymphocyte antigen 17

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Loaded Recombinant Mouse CD32b(C-6His) (RPES0629) on HIS1K Biosensor, can bind Recombinant Mouse IgG1 Fc (RPES1098) with an affinity constant of 0.388 uM as determined in BLI assay.

**Calculated MW:** 21.6 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** 33 kDa

**Reconstitution:** Please refer to the printed manual for detailed information. Loaded Recombinant Mouse

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

CD32b(C-6His) (RPES0629) on HIS1K Biosensor, can bind Recombinant Mouse IgG1 Fc (RPES1098) with an affinity constant of 0.388 uM as determined in BLI assay.

**Accession:** P08101

**Source:** HEK293 Cells-derived Mouse CD32b/FCGR2B protein Thr30-Pro210, with an C- terminal His

**Sequence:** Thr30-Pro210

**Form:** Lyophilized powder

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.