

# Recombinant Human ERK3/MAPK12 Protein (His & GST Tag)

RPES0691

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES0691

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** ERK-6, ERK3, ERK6, MAPK12, P38-gamma, P38GAMMA, PRKM12, SAPK-3, SAPK3

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His-GST

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** Baculovirus-Insect Cells

**Bio-Activity:** Not validated for activity

**Calculated MW:** 69.8 kDa

**Formulation:** Supplied as sterile solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH 8.0

**Observed MW:** 65 kDa

**Reconstitution:** -

**Accession:** P53778

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** Baculovirus-Insect Cells-derived Human ERK3/MAPK12 protein Met 1-Leu 367, with an N-terminal His & GST

**Sequence:** Met 1-Leu 367

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Form:** Liquid