

Recombinant Mouse FGF-9/FGF9 Protein (His Tag)

RPES0738

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0738

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Mouse

Endotoxin: < 0.01 EU per µg of the protein as determined by the LAL method.

Synonyms: FGF-9, Fgf-9, Fgf9, Fibroblast growth factor 9, GAF, Glia-activating factor, HBGF-9, heparin-binding growth factor-9

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 24.4 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 5% Trehalose, 1mM EDTA, 20% Glycerol, 1mM DTT, pH 8.5.

Observed MW: 25 kDa

Reconstitution: -

Accession: P54130

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Mouse FGF-9/FGF9 protein Met1-Ser208, with an N-terminal His

Sequence: Met1-Ser208

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid