

Recombinant Mouse PDGF-BB Protein (His Tag)

RPES0906

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0906

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Becaplermin, PDGF Subunit B, PDGF-2, PDGF2, PDGFB, Platelet-Derived Growth Factor B Chain, Platelet-Derived Growth Factor Beta Polypeptide, Platelet-Derived Growth Factor Subunit B, Proto-Oncogene c-Sis, SIS

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: E.coli

Bio-Activity: Measured in a cell proliferation assay using Balb/c3T3 mouse fibroblast cells. The ED 50 for this effect is 38.7 ng/ml.

Calculated MW: 13.4 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of 4mM HCl. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 15 kDa

Reconstitution: Please refer to the printed manual for detailed information. Measured in a cell proliferation assay using Balb/c3T3 mouse fibroblast cells. The ED50 for this effect is 38.7 ng/ml.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: AAH53430.1

Source: E.coli-derived Mouse PDGF-BB protein Ser82-Thr190, with an C-terminal His

Sequence: Ser82-Thr190

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.