

# Recombinant Human IA2/PTPRN Protein (aa 607-Asn686 &aa 795-888, RPES0916

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES0916

**Purity:** > 90 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** ICA 512, ICA3, ICA512, Islet cell antigen 512, Islet cell autoantigen 3, PTP IA-2, PTPRN, R-PTP-N, Receptor-type tyrosine-protein phosphatase-like N

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 22.9 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 1mM DTT, 1mM EDTA, pH 8.0.

**Observed MW:** 26-30 kDa

**Reconstitution:** -

**Accession:** Q16849

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human IA2, PTPRN protein Gln607-Asn686&Trp795-Leu888, with an N- terminal His

**Sequence:** Gln607-Asn686&Trp795-  
Leu888

**Form:** Liquid

**Notes:** Centrifuge before opening to ensure complete  
recovery of vial contents.