

Recombinant Mouse Granzyme B(G,H)/GZMB Protein (His Tag)

RPES0956

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES0956

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: CCP1, CTLA-1, Ctla-1, Ctla1, Cytotoxic cell protease 1, Fragmentin-2, Granzyme B(G,H), Gzmb

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 26.5 kDa

Formulation: Supplied as a 0.2 µm filtered solution of PBS, pH7.4.

Observed MW: 35 kDa

Reconstitution: -

Accession: P04187

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: HEK293 Cells-derived Mouse Granzyme B/GZMB protein Gly19-Ser247, with an C- terminal His

Sequence: Gly19-Ser247

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid