

# Recombinant Human PAFAH1B2/PAFAHB Protein (His Tag)

RPES1095

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1095

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** PAF Acetylhydrolase 30 kDa Subunit, PAF-AH 30 kDa Subunit, PAF-AH Subunit Beta, PAFAH Subunit Beta, PAFAH1B2, PAFAHB, Platelet-Activating Factor Acetylhydrolase IB Subunit Beta

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 26.6 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

**Observed MW:** 31 kDa

**Reconstitution:** -

**Accession:** P68402

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human PAFAH1B2, PAFAHB protein Ser2-Ala229, with an C-terminal His

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Ser2-Ala229

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid