

Recombinant Mouse IL1R2/CD121b Protein (Fc Tag)

RPES1193

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1193

Purity: > 80 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Interleukin-1 receptor type, CD121 antigen-like family member B, CD121b, IL-1 type II receptor, Il1r-2, IL-1R-2, IL-1R-beta, IL-1RT2, IL-1RT-2, Interleukin-1 receptor beta, Interleukin-1 receptor type 2, Interleukin-1 receptor type II, Il-1r2, Il1rb, Il1r2, IL-1R type II, CDw121b, CDw121b, IL1R2c, Interleukin 1 receptor, interleukin 1 receptor type 2, type II

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-Fc

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 64.9 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 85-110 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: P27931

Source: HEK293 Cells-derived Mouse IL1R2/CD121b protein Phe14-Glu355 , with an C-terminal Fc

Sequence: Phe14-Glu355

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.