

Recombinant Human TGFB3 Protein

RPES1241

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1241

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 0.01 EU per µg of the protein as determined by the LAL method.

Synonyms: LAP, Latency-associated peptide, TGF-beta-3, TGFB3, Transforming growth factor beta-3

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Measured by its ability to inhibit the IL-4-dependent proliferation of TF-1 mouse T cells. The ED50 for this effect is 10-80 pg/ml.

Calculated MW: 12.7 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of 50mM Glycine-HCl, 150mM NaCl, pH 2.5. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 12-14 kDa

Reconstitution: Please refer to the printed manual for detailed information. Measured by its ability to inhibit

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

the IL-4-dependent proliferation of TF-1 mouse T cells.
The ED50 for this effect is 10-80 pg/ml.

Accession: P10600

Source: HEK293 Cells-derived Human
TGFB3 protein Ala301-Ser412(Tyr340Phe)

Sequence: Ala301-Ser412(Tyr340Phe)

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.