

Recombinant Human SULT2A1 Protein (His Tag)

RPES1404

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1404

Purity: > 97 % as determined by reducing SDS-PAGE.

Contents: 50µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: Please contact us for more information.

Synonyms: Bile Salt Sulfotransferase, DHEA-ST, DHEAS, Dehydroepiandrosterone Sulfotransferase, HST, Hydroxysteroid Sulfotransferase, ST2, ST2A1, ST2A3, STD, SULT2A1, Sulfotransferase 2A1

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 34.5 kDa

Formulation: Lyophilized from sterile 50mM Tris, 500mM NaCl, 20% glycerol, pH 8.0 Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 35 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Accession: NP_003158.2

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in

Source: E.coli-derived Human SULT2A1 protein Ser 2-Glu 285, with an N-terminal His

Sequence: Ser 2-Glu 285

Form: Lyophilized powder

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.