

# Recombinant Human TLR4/CD284 Protein (His Tag)

RPES1421

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1421

**Purity:** > 87 % as determined by reducing SDS-PAGE.

**Contents:** 100µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** ARMD10, CD284, TLR-4, TLR4, TOLL

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** Baculovirus-Insect Cells

**Bio-Activity:** Immobilized Recombinant Human TLR4 / CD284 Protein (His Tag) (RPES1421) at 1 µg/mL (100 µL/well) on His Tag Antibody, Mouse MAb precoated (5 µg/mL, 100 µL/well) can bind Anti-TLR4 Antibody, the EC50 is 4-15 ng/mL.

**Calculated MW:** 70.5 kDa

**Formulation:** Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** 68 kDa

**Reconstitution:** Please refer to the printed manual for detailed information. Immobilized Recombinant

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

Human TLR4 / CD284 Protein (His Tag) (RPES1421) at 1 µg/mL (100 µL/well) on His Tag Antibody, Mouse MAb precoated (5 µg/mL, 100 µL/well) can bind Anti-TLR4 Antibody, the EC50 is 4-15 ng/mL. Immobilized Anti-TLR4 Antibody at 2 µg/mL (100 µL/well) can bind Recombinant Human TLR4 / CD284 Protein (His Tag) (RPES1421), the EC50 is 4-15 ng/mL .

**Accession:** O00206-1

**Source:** Baculovirus-Insect Cells-derived Human TLR4/CD284 protein Met 1-Lys631, with an C-terminal His

**Sequence:** Met 1-Lys631

**Form:** Lyophilized powder

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.