

# Recombinant Mouse IL1R1/CD121a Protein (His Tag)

RPES1533

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1533

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Mouse

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** D2S, Interleukin-1 receptor type, IL1-R, mL-1R, CD121 antigen-like family member A, CD121a, CD121b, D2S1473, IL-1 RI, IL1R, IL1r-1, IL1-R1, IL-1R-1, IL-1R-alpha, IL-1RT1, IL-1RT-1, IL-iR, Interleukin-1 receptor alpha, Interleukin-1 receptor type 1, mL-1R1, P80, IL-1r1, IL1ra, IL1r1, IL-1R type I, IL1R, Interleukin 1 receptor, interleukin 1 receptor type 1, mL-1R1, type I

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Not validated for activity

**Calculated MW:** 38.1 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** 50-90 kDa

**Reconstitution:** Please refer to the printed manual for detailed information.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Accession:** P13504

**Source:** HEK293 Cells-derived Mouse IL1R1/CD121a protein Leu20-Lys338, with an C-terminal His

**Sequence:** Leu20-Lys338

**Form:** Lyophilized powder

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.