

# Recombinant Human Cyclophilin D/PPID Protein (His Tag)

RPES1544

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1544

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** 40 kDa Peptidyl-Prolyl Cis-Trans Isomerase, CYP-40, CYP40, CYPD, Cyclophilin-40, Cyclophilin-Related Protein, PPID, PPIase D, Peptidyl-Prolyl Cis-Trans Isomerase D, Rotamase D

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His & C-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 43.9 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of PBS, 10% Glycerol, pH 7.4.

**Observed MW:** 40 kDa

**Reconstitution:** -

**Accession:** Q08752

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human Cyclophilin D, PPID protein Met 1-Ala370, with an N-terminal His & C-terminal His

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Met 1-Ala370

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid