

# Recombinant Human FKBP3/FKBP25 Protein (His Tag)

RPES1583

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1583

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** 25 kDa FK506-binding protein, 25 kDa FKBP, FK506-binding protein 3, FKBP-25, FKBP-3, FKBP25, Immunophilin FKBP25, PPIase FKBP3, Peptidyl-prolyl cis-trans isomerase FKBP3, Rapamycin-selective 25 kDa immunophilin, Rotamase

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 27.3 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 1mM DTT, 10% Glycerol, pH 8.0.

**Observed MW:** 30 kDa

**Reconstitution:** -

**Accession:** Q00688

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human FKBP3, FKBP25 protein Met 1-Asp224, with an N-terminal His

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Met 1-Asp224

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid