

Recombinant Human Cystathionine γ -Lyase/CTH Protein

RPES1619

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1619

Purity: > 85 % as determined by reducing SDS-PAGE.

Contents: 50 μ g, 10 μ g
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per μ g of the protein as determined by the LAL method.

Synonyms: CTH, Cystathionine Gamma-Lyase, Cysteine-Protein Sulphydrase, Gamma-Cystathionase

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 44.7 kDa

Formulation: Supplied as a 0.2 μ m filtered solution of 20mM Tris-HCl, 8% Sucrose, 0.05% Tween 80, pH8.0.

Observed MW: 38-50 kDa

Reconstitution: -

Accession: P32929

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human Cystathionine γ -Lyase, CTH protein Met 1-Ser405

Sequence: Met 1-Ser405

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid