

# Recombinant Human IL1R1/CD121a Protein (His Tag)

RPES1676

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1676

**Purity:** > 97 % as determined by reducing SDS-PAGE.

**Contents:** 100µg, 1mg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** D2S, Interleukin-1 receptor type, IL1RT, IL1-R, IL1R1, CD121A, D2S1473, IL-1R-alpha, IL1R, IL1RA, P80, CD121 antigen-like family member A, IL-1 RI, IL1-R1, IL-1R1, IL-1R-1, IL-1RT1, IL-1RT-1, Interleukin-1 receptor alpha, Interleukin-1 receptor type 1, IL1RT1, IL-1R type I, IL-1RI

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Measured by its ability to inhibit Interferon gamma secretion by human natural killer lymphoma NK-92 cells in the presence of 500 pg/mL IL1b. The EC50 for this effect is typically 40-200 ng/mL.

**Calculated MW:** 37.7 kDa

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** 55-57 kDa

**Reconstitution:** Please refer to the printed manual for detailed information.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Accession:** NP\_000868.1

**Source:** HEK293 Cells-derived Human IL1R1/CD121a protein Met 1-Thr 332, with an C-terminal His

**Sequence:** Met 1-Thr 332

**Form:** Lyophilized powder

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.