

Recombinant Human IL10-RA/IL-10 R α Protein (His Tag)

RPES1770

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1770

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50 μ g, 10 μ g
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU per μ g of the protein as determined by the LAL method.

Synonyms: IL-10R subunit, Interleukin-10 receptor subunit, IL-10R, IL10RA, CD210, CD210a, CDW210A, HIL-10R, IL-10R1, IL10R, IL-10 receptor subunit alpha, IL-10R subunit 1, IL-10R subunit alpha, IL-10RA, Interleukin-10 receptor subunit 1, Interleukin-10 receptor subunit alpha, IL-10R1, IL-10Ra, CDw210a antigen, HIL 10R, HIL10R, I10R1, IL 10R, IL 10R 1, IL 10R A, IL 10R1, IL10R1, Interleukin 10 receptor alpha, Interleukin 10 receptor alpha chain

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 25.2 kDa

Formulation: Lyophilized from a 0.2 μ m filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 38-59 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: Q13651

Source: HEK293 Cells-derived Human IL10-RA, IL-10 Ra protein His22-Asn235, with an C- terminal His

Sequence: His22-Asn235

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.