

Recombinant Mouse Interleukin-22/IL-22 Protein (mFc Tag)

RPES1785

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES1785

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Interleukin, Cytokine Zcyto, IL-D, TIFIL, Cytokine Zcyto18, IL-10-related T-cell-derived-inducible factor, IL-22a, IL-D110, IL-TIF, Interleukin-22a, TIFa, TIFIL-23, UNQ3099/PRO10096, zcyto18, Interleukin-22, T-cell-derived-inducible factor(IL-TIF), Interleukin-22a(IL-22a), IL-TIF alpha, IL-10-related, IL-22, Il22a, Iltif, Iltifa, Il22, TIFa

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-mFc

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 43.8 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 50-65 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: Q9JJY9

Source: HEK293 Cells-derived Mouse Interleukin-22/IL-22 protein Leu34-Val179, with an C- terminal mFc

Sequence: Leu34-Val179

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.