

# Recombinant Human MKK6 Protein (207 Ser/Asp, 211 Thr/Asp, His RPES1910

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES1910

**Purity:** > 88 % as determined by reducing SDS-PAGE.

**Contents:** 50µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** MAP2K6 , MAPKK6, MEK6, MKK6, PRKMK6, SAPKK-3, SAPKK3

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His-GST

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** Baculovirus-Insect Cells

**Bio-Activity:** The specific activity was determined to be 1250 nmol/min/mg using inactive MAPK14 as substrate.

**Calculated MW:** 65.3 kDa

**Formulation:** Supplied as sterile solution of 20mM Tris, 500mM NaCl, pH 8.0, 10% glycerol

**Observed MW:** 60 kDa

**Reconstitution:** -

**Accession:** P52564-1

**Source:** Baculovirus-Insect Cells-derived Human MKK6 protein Met 1-Asp 334, Ser 207/Asp, Thr 211/Asp, with an N-terminal His & GST

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Met 1-Asp 334; Ser 207/Asp;  
Thr 211/Asp  
**Form:** Liquid

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.