

# Recombinant Human UBE2T Protein (His Tag)

RPES2090

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES2090

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** Cell Proliferation-Inducing Gene 50 Protein, UBE2T, Ubiquitin Carrier Protein T, Ubiquitin-Conjugating Enzyme E2 T, Ubiquitin-Protein Ligase T

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 24.7 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

**Observed MW:** 25-30 kDa

**Reconstitution:** -

**Accession:** Q9NPD8

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human UBE2T protein Met 1-Val197, with an N-terminal His

**Sequence:** Met 1-Val197

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid