

Recombinant Human UBE2Z Protein (His Tag)

RPES2190

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2190

Purity: > 90 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: HOYS7, UBE2Z, Uba6-Specific E2 Conjugating Enzyme 1, Ubiquitin Carrier Protein Z, Ubiquitin-Conjugating Enzyme E2 Z, Ubiquitin-Protein Ligase Z, Use1

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 30.2 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 50mM HEPES, 50mM NaCl, pH 8.0.

Observed MW: 32 kDa

Reconstitution: -

Accession: Q9H832-2

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human UBE2Z protein Met 1-Val246, with an N-terminal His

Sequence: Met 1-Val246

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid