

Recombinant Human AKR1C2 Protein

RPES2268

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2268

Purity: > 90 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: 2-Dihydrobenzene-1, 2-Diol Dehydrogenase, 3-Alpha-HSD3, Aldo-Keto Reductase Family 1 Member C2, Chlordecone Reductase Homolog HAKRD, DD-2, DD/BABP, DD2, Dihydrodiol Dehydrogenase 2, Dihydrodiol Dehydrogenase/Bile Acid-Binding Protein, Trans-1, Type III 3-Alpha-Hydro

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 36.7 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, pH 8.0.

Observed MW: 35 kDa

Reconstitution: -

Accession: P52895

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human AKR1C2 protein Met 1-Tyr323

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Met 1-Tyr323

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid