

Recombinant Human GALNT3 Protein (His Tag)

RPES2332

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2332

Purity: > 85 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: GalNAc-T3, HFTC, HHS, Polypeptide GalNAc transferase 3, Polypeptide N-acetylgalactosaminyltransferase 3, Protein-UDP acetylgalactosaminyltransferase 3, UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 3, pp-GaNTase 3

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 69.1 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Observed MW: 80 kDa

Reconstitution: -

Accession: Q14435

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in

Source: HEK293 Cells-derived Human GALNT3 protein Gln38-Asp633, with an C-terminal His

Sequence: Gln38-Asp633

Form: Liquid

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.