

Recombinant Human PPM1G/PP2C-gamma Protein (aa 317-546, His Tag)

RPES2411

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2411

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 100µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: Please contact us for more information.

Synonyms: PP2C-Gamma, PP2CG, PP2CGAMMA, PPM1G, Protein Phosphatase 1C, Protein Phosphatase 1G, Protein Phosphatase 2C Isoform Gamma, Protein Phosphatase Magnesium-Dependent 1 Gamma

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 26.6 kDa

Formulation: Lyophilized from sterile 50mM Tris, 1mM DTT, 20% glycerol, pH 7.5 Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 30 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Accession: O15355

Source: E.coli-derived Human PPM1G/PP2C-gamma protein Met 317-Asp 546, with an N- terminal His

Sequence: Met 317-Asp 546

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.