

Recombinant Human GRK6/GPRK6 Protein (His & GST Tag)

RPES2620

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2620

Purity: > 93 % as determined by reducing SDS-PAGE.

Contents: 50µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: GPRK6

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His-GST

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: Baculovirus-Insect Cells

Bio-Activity: The specific activity was determined to be 7 nmol/min/mg using casein as substrate.

Calculated MW: 95.1 kDa

Formulation: Supplied as sterile solution of 20mM Tris, 500mM NaCl, 2mM GSH, 0.5mM PMSF, pH 7.4

Observed MW: 85 kDa

Reconstitution: -

Accession: P43250-2

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: Baculovirus-Insect Cells-derived Human GRK6/GPRK6 protein Met 1-Arg 589, with an N-terminal His & GST

Sequence: Met 1-Arg 589

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid