

Recombinant Human KLK-8/Kallikrein-8 Protein (His Tag)

RPES2890

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES2890

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: KLK8, Kallikrein-8, NP, NRPN, Neuropsin, Ovasin, PRSS19, Serine Protease 19, Serine Protease TADG-14, TADG14, Tumor-Associated Differentially Expressed Gene 14 Protein, hK8

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 26.0 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.

Observed MW: 35 kDa

Reconstitution: -

Accession: AAH40887.1

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: HEK293 Cells-derived Human KLK-8, Kallikrein-8 protein Gln29-Gly260, with an C- terminal His

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Gln29-GLy260

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid