

# Recombinant Human SH2D1A/SAP Protein (His Tag)

RPES3073

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES3073

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** DSHP, Duncan Disease SH2-Protein, EBVS, IMD5, LYP, MTCP1, SAP/SH2D1A, SH2 Domain-Containing Protein 1A, SH2D1A, SLAM-Associated Protein, Signaling Lymphocytic Activation Molecule-Associated Protein, T-Cell Signal Transduction Molecule SAP, XLP, XLPD, XLPD1

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 16.4 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 7.5.

**Observed MW:** 16 kDa

**Reconstitution:** -

**Accession:** O60880

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Source:** E.coli-derived Human  
SH2D1A/SAP protein Met 1-Pro128, with  
an N-terminal His

**Sequence:** Met 1-Pro128

**Form:** Liquid

**Notes:** Centrifuge before opening to ensure complete  
recovery of vial contents.