

Recombinant Human ATG10 Protein (His&T7 Tag)

RPES3103

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES3103

Purity: > 90 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: APG10-Like, ATG10, Autophagy-Related Protein 10, Ubiquitin-Like-Conjugating Enzyme ATG10

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-T7 & C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 24.5 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 8% Sucrose, 20% Glycerol, 0.02% Tween 80, 3mM TCEP, pH 8.5.

Observed MW: 30 kDa

Reconstitution: -

Accession: Q9H0Y0

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human ATG10 protein Met 1-Thr190, with an N-terminal T7 & C- terminal His

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Met 1-Thr190

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid