

Recombinant Human MAN1B1 Protein (His Tag)

RPES3201

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES3201

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: 2-Alpha-Mannosidase, 2-Mannosidase, ER Alpha-1, ER Mannosidase 1, ERMan1, Endoplasmic Reticulum Mannosyl-Oligosaccharide 1, MAN1B1, Man9GlcNAc2-Specific-Processing Alpha-Mannosidase, Mannosidase Alpha Class 1B Member 1

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 68.7 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 50mM Tris-HCl, 10mM reduced Glutathione, pH 8.0.

Observed MW: 58-80 kDa

Reconstitution: -

Accession: Q9UKM7

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: HEK293 Cells-derived Human MAN1B1 protein Asp106-Ala699, with an C-terminal His

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Asp106-Ala699

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid