

Recombinant Human LEFTY2 Protein (His Tag)

RPES3222

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES3222

Purity: > 90 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: EBAF, Endometrial bleeding-associated factor, LEFTA, LEFTY2, LEFTYA, Left-right determination factor 2, Left-right determination factor A, Protein lefty-2, Protein lefty-A, TGF-beta-4, TGFB4, Transforming growth factor beta-4

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 39.1 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of 20mM Histidine-HCl, 4% Sucrose, 4% Mannitol, 0.1% Tween 80, pH6.0. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 45-50 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Accession: O00292

Source: HEK293 Cells-derived Human LEFTY2 protein Phe78-Pro366, with an N-terminal His

Sequence: Phe78-Pro366

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.