

# Recombinant Human IL17RA Protein (Fc Tag)

RPES3390

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES3390

**Purity:** > 90 % as determined by reducing SDS-PAGE.

**Contents:** 50µg, 10µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** IMD, MGC, CTLA, CANDF, CDw, IL17RA, CANDF5, CD217, CDw217, IL-17RA, IL17R, IMD51, hIL-17R, CD127, CD217 antigen, CDW127, CTLA8, HIL 17R, hIL17R, I17RA, IL 17 receptor, IL 17 receptor A, IL 17RA, IL17, IL-17 R alpha, IL-17 receptor A, IL17A receptor, IL-17R, IL17R-A, IL17-RA, IL-17R-alpha, IL7RA, ILRA, Interleukin 17 receptor A, Interleukin-17 receptor A, MGC10262, CDw217, Interleukin-17 receptor A

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-Fc

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Immobilized Human IL-17RA-Fc at 2µg/ml (100 µl/well) can bind Anti-Human IL-17RA mAb. The ED50 of Anti-Human IL-17RA mAb is 13.50 ng/ml.

**Calculated MW:** 60.6 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Observed MW:** 90-120 kDa

**Reconstitution:** Please refer to the printed manual for detailed information.

**Accession:** Q96F46

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** HEK293 Cells-derived Human IL17RA protein Leu33-Trp320, with an C-terminal Fc

**Sequence:** Leu33-Trp320

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Lyophilized powder