

Recombinant Mouse TGFB2 Protein

RPES3426

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES3426

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Mouse

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: BSC-1 cell growth inhibitor, Cetermin, G-TSF, Glioblastoma-derived T-cell suppressor factor, MGC116892, Polyergin, TGF-beta-2, TGF-beta2, TGFB2, transforming growth factor beta-2

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Immobilized Mouse TGFB2 at 10µg/ml(100 µl/well) can bind Human TGFBR2-FC(RPES4049). The ED50 of Mouse TGFB2 is 0.136µg/mL.

Calculated MW: 12.7 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of 4mM HCl. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 12 kDa

Reconstitution: Please refer to the printed manual for detailed information. Immobilized Mouse TGFB2 at 10µg/ml(100 µl/well) can bind Human TGFBR2-

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

FC(RPES4049). The ED50 of Mouse TGFB2 is 0.136µg/mL.

Accession: P27090

Source: HEK293 Cells-derived Mouse TGFB2 protein Ala303-Ser414

Sequence: Ala303-Ser414

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.