

Recombinant E.coli Beta-galactosidase Protein

RPES3494

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES3494

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: E.coli

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Beta-gal, Beta-galactosidase, Lactase, lacZ

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 112.9 kDa

Formulation: Supplied as a 0.2 µm filtered solution of PBS, pH7.4.

Observed MW: 115 kDa

Reconstitution: -

Accession: P00722

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived E.coli Beta-galactosidase protein Met1-Lys1024(12-41AA deletion)

Sequence: Met1-Lys1024(12-41AA deletion)

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid