

# Recombinant Human GFRA1/GDNFRA Protein (aa 1-424, His Tag)

RPES4321

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES4321

**Purity:** > 98 % as determined by reducing SDS-PAGE.

**Contents:** 100µg, 1mg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** GDNF Family Receptor Alpha-1, GDNF Receptor Alpha-1, GDNFR-Alpha-1, GDNFRA, GFR-Alpha-1, GFRA1, RET Ligand 1, RETL1, TGF-Beta-Related Neurotrophic Factor Receptor 1, TRNR1

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** 1. Measured in a cell proliferation assay using SH- SY5Y human neuroblastoma cells. The ED50 for this effect is typically 0.2-1 µg/mL in the presence of 40 ng/mL Recombinant Human GDNF. 2. Measured by its ability to bind human GDNF (native) in functional ELISA.

**Calculated MW:** 46.0 kDa

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Observed MW:** 55-60 kDa

**Reconstitution:** Please refer to the printed manual for detailed information.

**Accession:** NP\_665736.1

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** HEK293 Cells-derived Human GFRA1/GDNFRA protein Met 1-Ser 424, with an C- terminal His

**Sequence:** Met 1-Ser 424

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Lyophilized powder