

Recombinant Human BCL2L1/Bcl-XL Protein (His Tag)

RPES4389

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES4389

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Apoptosis Regulator Bcl-X, BCL-XL/S, BCL2L, BCLX, BCLXL, BCLXS, Bcl-2-Like Protein 1, Bcl-X, Bcl2-L-1, PPP1R52, bcl-xL, bcl-xS

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 24.9 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM HEPES, 50mM KCl, 20% Glycerol, pH 7.5.

Observed MW: 29 kDa

Reconstitution: -

Accession: Q07817

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human BCL2L1, Bcl-XL protein Met1-Arg212, with an C-terminal His

Sequence: Met1-Arg212

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Form: Liquid