

Recombinant Human LAMTOR3/MP1 Protein (His Tag)

RPES4401

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES4401

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: LAMTOR3, Late Endosomal/Lysosomal Adaptor and MAPK and MTOR Activator 3, MAP2K1IP, MEK-Binding Partner 1, Mitogen-Activated Protein Kinase Kinase 1-Interacting Protein 1, Mitogen-Activated Protein Kinase Scaffold Protein 1, Mp1, Ragulator Complex Protein LAMTOR3

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 15.8 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, 5mM DTT, 30% Glycerol, pH 8.0.

Observed MW: 15 kDa

Reconstitution: -

Accession: Q9UHA4

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Source: E.coli-derived Human LAMTOR3, MP1 protein Met 1-Val124, with an N-terminal His

Sequence: Met 1-Val124

Form: Liquid

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.