

Recombinant Human Alpha-Synuclein/SNCA Protein

RPES4685

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES4685

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: Alpha-Synuclein, NACP, Non-A Beta Component of AD Amyloid, Non-A4 Component of Amyloid Precursor, PARK1, SNCA

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: None

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: E.coli

Bio-Activity: Immobilized Recombinant Human SNCA(RPES4685) at 2µg/ml (100 µl/well) can bind Anti-Human SNCA antibody. The ED50 of Anti-SNCA antibody is 3.11ng/ml.

Calculated MW: 14.5 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 17 kDa

Reconstitution: Please refer to the printed manual for detailed information. Immobilized Recombinant Human SNCA(RPES4685) at 2µg/ml (100 µl/well) can

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

bind Anti-Human SNCA antibody. The ED50 of Anti-SNCA antibody is 3.11 ng/ml.

Accession: P37840

Source: E.coli-derived Human SNCA protein Met1-Ala140

Sequence: Met1-Ala140

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.